


Leading the way in microfluidics in Point-of-care solutions

Diagnostics is a dynamic, growing industry that relies on faster and more accurate results to improve clinical outcomes, increase operational efficiencies, and reduce overall costs.





You are innovating


- To diagnose and detect nearer, faster, and cheaper
- To reduce manual handling steps and complexity
- To mitigate misdiagnosis

To ideally create

A device that is accurate, sensitive, rapid, and affordable to advance medical diagnosis¹

Micro functions make all the difference in microfluidics

Studies in microfluidics find that air vents accomplish a number of versatile functions than just air venting. Yet, as a standalone component, enables portability for the point-of-care device.



vent & filter

Allows air expression and pressurization; control of fluid delivery and metering; mitigate bubbling; eliminates sample evaporation; minimizes humidity; while acting as a bacterial and viral barrier to filter out contamination

Porex Virtek®
hydrophobic PTFE vents
make all the difference in
microfluidics

Agenda

1. Support fluid flow and pressure equalization
2. Maintain a liquid tight, closed system barrier to minimize leakages and biohazards
3. Maintain sterility and minimize outside contaminants for sample integrity
4. A lower MVTR material to protect dried enzyme-based reagents from early reconstruction and minimize sample evaporation for test accuracy
5. A proven and effective mechanism for bubbling challenges

A Porex Virtek® Hydrophobic PTFE vent for fluid control

Keep fluid inside

Omni-directional pore structure with no front or back orientations, so membrane functions to **keep fluids in**, yet allows air to escape.

Air displacement

Air displacement supports **fluid movement**, which allows it to align and be accurately metered in volume in some designs¹.

Pressure equalization


Proper air venting management, for **pressure equalization** during heating and cooling cycles.





A Porex Virtek® Hydrophobic PTFE vent to maintain a closed system

- **Prevents leakage** and exposure of reagent amplicons where user handling errors may occur.
- Seals liquids inside, to **manage biological waste** after disposal to minimize spread of biohazards, especially with infectious disease testing.
- Surface **readily vents after** contact with liquid, unlike other porous membrane materials which can become blocked after fluid contact.
- Both sides **perform the same** unlike many cast or stretched membranes that have non-functional supporting layers.





Active venting maintains sensitivity by preventing contaminants entering the device, critical to sample integrity

Porex Virtek® PTFE has extremely high VFE & BFE (up to log 6) for high barrier protection from a durable vent with optimal airflow to meet (or exceed) FDA & EU guidelines.

A barrier against viral and bacterial particles¹

Bacterial filtration efficiency (BFE)²:

- Determines the filtration efficiency by comparing bacterial counts to test article effluent counts.
- BFE greater than 99.99% filtration efficiency demonstrates the material is effective in blocking bacterial micro-organisms.


Viral filtration efficiency (VFE)²:


- VFE testing takes BFE a step further by measuring virus-containing aerosol particles of specific sizes to determine efficiency of filtration media in capturing those aerosols.
- VFE greater than 99.99% filtration efficiency demonstrates the material is effective in preventing infectious pathogens entering or exiting from test device.

Up to **5 x lower MVTR** than stretched porous membranes for a liquid-tight barrier that is **critical to stability and sensitivity of test**¹

The **biggest roadblock** to commercialization and regulatory approval is accuracy, sensitivity, and specificity of the test². Porex Virtek® PTFE vent membranes **have up to 5 times lower Moisture Vapor Transmission Rates** than stretched porous membranes to:

- Provide a **liquid tight barrier** to protect lyophilized reagents and minimize moisture and humidity exposure, avoiding early reconstruction.
- **Preserve liquid reagents** stored within a sterile, vapor-resistant reaction chamber, also critical to sensitivity of test.
- **Minimizes sample evaporation** during temperature cycling.
- Supportive to **storage and shelf life** of enzyme-based reagents.





Porex Virtek[®] PTFE as a **passive "debubbling" mechanism** to mitigate bubbles that become trapped

- The unique membrane surface structure promotes **microbubbles to burst** and exit.
- Allows **air pass through** of the vent to promote fluid movement and transport bubbles out of the system.
- Keeping sample liquids in, **minimizing sample evaporation** and expulsion of reagents.
- Yet **does not block** at contact with liquid, unlike some other porous membrane materials.

Micro functions that mitigate the challenge of bubbling in microfluidics

Bubbles occur due to

- Portability of device when tilting and shaking during transportation¹
- Oxidation of surface tension and irregularities when bringing liquid inside of microchannels²
- During thermal cycling or changes in pressure conditions³
- Commonly used chip plate materials like PDMS acrylic copolymer have relatively low gas removal which influences bubble formation⁴

Creating issues such as

- Difficult to remove micro-bubbles trapped within reaction chambers, metering channels, and optical systems, impacting quantitative measurements⁵
- Increases in formation and size with changes in temperature and pressure conditions⁶
- Distortion of fluid flow⁷
- Causes damage to cells at liquid-gas interface⁷
- Causes evaporation and expulsion of PCR reagents⁸

1. Combes RD, Balls M, Bhogal N. New technology for toxicity testing. *Adv Exp Med Biol*. 2012;745:v-xiii, xv, xvii passim.

2. Cheng, Hao-Bin & Lu, Yen-Wen. (2014). Applications of textured surfaces on bubble trapping and degassing for microfluidic devices. *Microfluidics and Nanofluidics*.

3. Prakash M, Gershenfeld N. Microfluidic bubble logic. *Science*. 2007 Feb 9;315(5813):832-5.

4. Sung JH, Shuler ML. Prevention of air bubble formation in a microfluidic perfusion cell culture system using a microscale bubble trap. *Biomed Microdevices*. 2009 Aug;11(4):731-8.

5. Podczerviensk, McDowell, Levine. Affect of Air Bubbles on Filling and Metering in a Microfluidic Device. *NSTI-Nanotech 2012*. Vol 2, 2012.


6. Liu, H. B., Gong, H. Q., Ramalingam, N., Jiang, Y., Dai, C. C., & Hui, K. M. (2007). Micro air bubble formation and its control during polymerase chain reaction (PCR) in polydimethylsiloxane (PDMS) microreactors. *Journal of Micromechanics and Microengineering*, 17(10), 2055–2064.

7. Liu, C., Thompson, J. and Bau, H. (2011) "A membrane-based, high-efficiency, microfluidic debubbler", *Lab on a Chip*, 11(9), p. 1688.

8. Trung, N. B., Saito, M., Takabayashi, H., Viet, P. H., Tamiya, E., & Takamura, Y. (2010). Multi-chamber PCR chip with simple liquid introduction utilizing the gas permeability of polydimethylsiloxane. *Sensors and Actuators B: Chemical*, 149(1), 284–290.

Porex Virtek® sintered PTFE is compatible for microfluidic applications

- **Naturally hydrophobic**
 - 100% pure PTFE with no additives or treatments, that naturally resists water
 - No water drawn through membrane under typical vacuum and pressure conditions
- **Thermal stability**
 - Highly heat resistant up to 260° C
- **Chemical inertness**
 - Compatible with reagents that often contain alcohols, surfactants, proteins, and salts
- **Robust, pure and durable**
 - Requires no supporting layers or chemical treatments
 - Safe to handle without damaging
 - No assembly orientation required
 - High purity with virtually no leachable or extractables
- **Manufactured in class 100k cleanroom**
- **Raw material certified**
 - USP class VI (bio-compatibility)
 - Free of PFOA



Cross section view




Free of PFOA

Compliant with (EC) 1907/2006 REACH / Regulation (EU) 2019/1021 POP


Microscopic comparison with other materials

Sintered PTFE Membrane




- Material has depth, is self supporting and will rebound under pressure
- Individual particles bonded to its neighbor for superior strength
- No additives or binders necessary
- Membrane is a depth filter and follows a torturous path
- Both sides identical and omnidirectional

Expanded / Stretched PTFE Membrane



- Active membrane has minimal depth and is only a surface filter
- Bonding to substrate requires additional adhesives or lamination to tie layers together
- Service temperatures limited to substrate
- Both sides have different properties
- Membrane properties can be altered by applied pressure

Cast Membrane




- Casting process often requires chemicals or other additives to create structure, leading to residuals and contamination issues
- Layers can delaminate
- Service temperatures and chemical resistance limited to substrate
- Membrane orientation is critical (different sides have different properties)
- Treatment often needed for hydrophobicity

Common **Assembly** options for Porex Virtek PTFE

There is **no standardization in design** for microfluidic devices, so sourcing and fitting custom components can be tricky. Porex Virtek® PTFE vents **make your process simpler by:**

- Easily heat sealed, vibrationally welded or numerous other assembly options
- Has no right-side assembly orientation
- Is easy to handle and does not easily damage or changes with contact
- Enables quicker production and device finishing in high-speed assembly



Material Property Ranges: physical properties

Material	Thickness mm nominal	Typical Airflow l/hr/cm ² at 70 mbar	Dry Filtration Efficiency* >99.99%	BFE** % Nominal	VFE*** % Nominal	Typical WEP**** mBar	WIP
MD10	0.13	125 (min 70)	0.5 µm	>99.99	~	270 (min 175)	350
MD10L	0.3	85 (min 48)	0.5 µm	>99.99	~	270 (min 175)	350
MD15	0.18	70 (min 45)	0.4 µm	>99.99	~	380 (min 265)	450
MD20	0.25	34 (min 16)	0.1 µm	>99.9999	>99.999	520 (min 350)	600
MD22	0.1	15 (min 5)	0.2 µm	>99.99	~	750 (min 500)	900
MD25	0.19	5 (min 2)	0.1 µm	>99.9999	~	1000 (min 750)	950

Looking for Pore Size? Unlike typical surface filter membranes, depth filtration membranes like Porex Virtek retain particles smaller than its nominal pore size through pore-size gradients and tortuous path. Use dry filtration efficiency to better compare.

Oleophobic treated materials available to repel low surface energy fluids (oils, alcohols & surfactants)

* According to IEST RP-CC007.2 2009

** Bacterial Filtration Efficiency (BFE) data is based on a modified version of ASTM F2101

*** Viral Filtration Efficiency (VFE), † Not tested but similar results to MD20 expected.

**** WEP (Water Entry Pressure)

Properties are typical and not meant for specifications. Selected options and adhesives may affect properties. Complete testing data and information is available upon request – all exceeded the standard BFE value of 98% RoHS, WEEE, REACH Compliant (PFOA Free)

Range of material options



- **Master roll:**
 - 330mm wide
 - Thickness from 0.1mm - 3.0mm
 - Roll length dependent on thickness (3m – 100m)
- **Converting options:**
 - **Slit Rolls**
 - 8mm minimum width (for most materials)
 - Provided on 76mm (3") ID plastic cores
 - **Die-cutting**
 - Minimum 3mm diameter
 - Minimum 5mm diameter if with adhesive
 - Custom Shapes
 - Robotic Frit cutting for thicker membranes
 - **Adhesive Discs**
 - Numerous stock, standard and custom size/shapes available
 - **Lamination**
 - PP / PE scrim options available
 - Adhesive backing



TOGETHER, WE ARE MAKING
THE WORLD SAFER, HEALTHIER
AND MORE PRODUCTIVE.