



## **Secondary Operations Part 1 - Solvent Bonding:** A versatile approach for fabricating of medical assemblies

Pooja Bajaj, PhD, Product Development Manager - Molding Compounds

Kay Bernhard, PhD, Technical Center Director - Molding Compounds

Andrew Sneeringer, Technical Marketing Specialist - Medical Molding Compounds

Maurice Biagini, Business Director - Medical Molding Compounds

# Agenda

- Introduction to CYROLITE®
- Why Solvent Bonding?
- Key factors influencing the bond strength – solvents, tubing materials, polymer etc.
- Mechanical testing to measure bond strength
- Overall performance of acrylic copolymers and other medical plastics

## CYROLITE® acrylic-based copolymer

# CYROLITE®

ADVANCED  
MEDICAL ACRYLICS

- ✓ **Leading acrylic copolymer** specifically developed for use in the Medical Device and Diagnostics Industries
- ✓ **Extensive development and application in healthcare** and labs for more than **40 years**
- ✓ **Offering eight grades** with balanced properties



# Global leader in medical grade acrylic copolymers

## Business Portfolio

### Infusion Therapy




- IV connectors
- Valves
- Y Sites
- Luer Locks
- IV Filter
- Catheter accessories

### Diagnostics



- Diagnostic cuvettes
- Rotors
- Test Packs
- Optical sensor view ports
- Trays etc.

### Blood Management



- Separators
- Filters
- Valves
- Connectors

### Medical Others



- Sterile packaging
- Respiratory components
- Ophthalmic components

Solvent bonding as a versatile fabrication technique  
for joining various plastic materials to medical tubing

# Promoting adhesion - Key fabrication strategies and factors

## Strategies to promote joining between polymers



Chemical bonding by solvents or epoxy




Pressure or Temperature sensitive



Eliminating adhesives via screws, inserts, or rivets




Surface modification by corona, plasma, UV, flame



Welding by ultrasonic or laser methods

## What drives Adhesion?



# Solvent Bonding

- **Motivation**

- Solvent Bonding is an established method to connect **CYROLITE®** to soft tubing using **compatible solvents**
- Incumbent **tubing material** of choice is DEHP PVC
- Industry is seeking non-DEHP PVC and PVC-free materials due to phthalate and incineration concerns (release of dioxins)

- **Objective**

- Improve bonding to non-DEHP PVC, PVC alternative solutions at  $\geq 11$  lb.f
- cost effective
- robust alternatives
- compliant for medical applications


- **Applications**

- Delivery systems, interconnectors, microfluidic devices




# Solvent Bonding

## Fabricating polymer - tubing assembly






# Choice of polymer and tubing material – Key considerations



*\*A review of alternatives to di (2-ethylhexyl) phthalate-containing medical devices in the neonatal intensive care unit. EDS Van Vliet, EM Reitano, JS Chhabra, GP Bergen, and RM Whyatt, J Perinatol. 2011, 31(8): 551–560.*


# Solvent application methods

## Fabricating polymer - tubing assembly




### Variables

- Solvent choice
- Solvent volume



Simple  
Cost - effective  
Strength of bond depends on compatibility




ESCR in amorphous polymers  
Extended drying times depending of solvent volatility

# How does a solvent influence bonding?

## Polymer compatibility - Hansen solubility parameter

Solubility parameters derived from the cohesive energies of polymer and solvent based on:

- $\delta_d$ , dispersion bonds
  - $\delta_p$ , polar bonds
  - $\delta_h$ , hydrogen bonds
- $R_a$  = distance between solvent and polymer
  - The smaller the distance  $R_a$  between solvent and polymer, the more the solvent will be able to dissolve the polymer
  - Relative energy difference  $RED = R_a/R_0$



*Schematic illustration of the Hansen space for a polymer and solvent compatibility*

# Solvent compatibility

## \*Solvents

<sup>1</sup>Cyclohexanone (CH)

<sup>2</sup>Methylene Chloride (MC)

Tetrahydrofuran


Toluene

Ethyl Acetate

90:10 CH (Cyclohexanone) / Methyl Ethyl Ketone (MEK)

Methyl Ethyl Ketone

Acetone




*\*Degree of overlap of spheres in Hansen space correlates with solvent - polymer compatibility*

1. 100% CH is reported to cause stress cracking in medical plastics
2. MC is reported to be carcinogenic, mutagenic and banned in many countries

# CYROLITE® GS-90

## Bond strength to non-DEHP PVC using compatible solvents



### \*Hansen Solubility Parameter $\delta$ ((MJ/m<sup>3</sup>)<sup>1/2</sup>)

| Substrate Polymer | $\delta D$ | $\delta P$ | $\delta H$ |
|-------------------|------------|------------|------------|
| Acrylics          | 18.1       | 8.4        | 3.1        |

### \*Hansen Solubility Parameter $\delta$ ((MJ/m<sup>3</sup>)<sup>1/2</sup>)

| **Solvents          | $\delta D$ | $\delta P$ | $\delta H$ |
|---------------------|------------|------------|------------|
| Cyclohexanone       | 17.8       | 8.4        | 5.1        |
| Methyl ethyl ketone | 15.9       | 8.7        | 5.1        |
| Acetone             | 15.5       | 10.4       | 6.9        |

\*Dispersion forces,  $\delta D$  ((MJ/m<sup>3</sup>)<sup>1/2</sup>)

Dipolar intermolecular force,  $\delta P$  ((MJ/m<sup>3</sup>)<sup>1/2</sup>)

Hydrogen bonds,  $\delta H$  ((MJ/m<sup>3</sup>)<sup>1/2</sup>)



\*Handbook of Solvents, Volume 1 - Properties (2nd Edition)

\*\*Dwivedi AM. Residual solvent analysis in pharmaceuticals. Pharmaceutical Technology. November (2002) 42


# Solvent bonding

## Fabricating CYROLITE® GS-90 - tubing assembly

Tubing to CYROLITE® Y-site connection




Joining assemblies




# Solvent bonding

## Fabricating CYROLITE® GS-90 - tubing assembly


Tubing to CYROLITE® GS-90 Y-site connection



Mechanical tests - Measure of bond strength



## Solvent bonding – Experimental parameters



### \*Solvents


Current

90:10 CH (Cyclohexanone) / Methyl Ethyl Ketone (MEK)




# CYROLITE® GS-90 shows strong bond strength to self and non-DEHP PVC to using Cyclohexanone/MEK (90/10)

## CYROLITE® GS-90 to CYROLITE® GS-90




## CYROLITE® GS-90 to non-DEHP PVC




Lap shear pull strength at 20inch/min (500mm/min)

# Thermoplastic Polyolefin tubing to CYROLITE® GS-90 shows weak bond strength using Cyclohexanone/MEK (90/10)




Lap shear pull strength at 20inch/min  
(500mm/min)



- Similar behavior is experienced with PC and MABS

# Thermoplastic Polyurethane (TPU) tubing to CYROLITE® GS-90 shows excellent bond strength using Cyclohexanone/MEK (90/10)



- Similar behavior is experienced with PC and MABS

## \*Hansen Solubility Parameter $\delta$ ((MJ/m<sup>3</sup>)<sup>1/2</sup>)

| Substrate Polymer | $\delta_D$ | $\delta_P$ | $\delta_H$ |
|-------------------|------------|------------|------------|
| Acrylics          | 18.1       | 8.4        | 3.1        |

## \*Hansen Solubility Parameter $\delta$ ((MJ/m<sup>3</sup>)<sup>1/2</sup>)


| Tubing material | $\delta_D$ | $\delta_P$ | $\delta_H$ |
|-----------------|------------|------------|------------|
| PVC             | 18.7       | 10         | 3.1        |
| TPU             | 18.8       | 10         | 8.2        |

## \*Hansen Solubility Parameter $\delta$ ((MJ/m<sup>3</sup>)<sup>1/2</sup>)

| Solvents            | $\delta_D$ | $\delta_P$ | $\delta_H$ |
|---------------------|------------|------------|------------|
| Cyclohexanone       | 17.8       | 8.4        | 5.1        |
| Methyl ethyl ketone | 15.9       | 8.7        | 5.1        |

# CYROLITE® **GS-90** bond superior to non-DEHP PVC compared to other medical plastics using Cyclohexanone/MEK (90/10)

## Bond strength (lbf)



## Hansen solubility parameter

| *Hansen Solubility Parameter $\delta$ ((MJ/m <sup>3</sup> ) <sup>1/2</sup> ) |            |            |            |
|------------------------------------------------------------------------------|------------|------------|------------|
| Substrate Polymer                                                            | $\delta_D$ | $\delta_P$ | $\delta_H$ |
| Acrylics                                                                     | 18.1       | 8.4        | 3.1        |
| MABS                                                                         | 18.5       | 3.7        | 3.9        |
| PC                                                                           | 19.6       | 8.8        | 5.7        |

| *Hansen Solubility Parameter $\delta$ ((MJ/m <sup>3</sup> ) <sup>1/2</sup> ) |            |            |            |
|------------------------------------------------------------------------------|------------|------------|------------|
| Tubing material                                                              | $\delta_D$ | $\delta_P$ | $\delta_H$ |
| PVC                                                                          | 18.7       | 10         | 3.1        |

| *Hansen Solubility Parameter $\delta$ ((MJ/m <sup>3</sup> ) <sup>1/2</sup> ) |            |            |            |
|------------------------------------------------------------------------------|------------|------------|------------|
| Solvents                                                                     | $\delta_D$ | $\delta_P$ | $\delta_H$ |
| Cyclohexanone                                                                | 17.8       | 8.4        | 5.1        |
| Methyl ethyl ketone                                                          | 15.9       | 8.7        | 5.1        |

\*Title: Handbook of Polymers (2nd Edition)  
Handbook of Solvents, Volume 1 - Properties (2nd Edition)

# CYROLITE® GS-90 and other Medical Plastics

## Solvent Bonding to non-DEHP PVC medical tubing and alternatives


| Solvent compatibility to non-DEHP PVC medical tubing | CYROLITE® GS-90 |
|------------------------------------------------------|-----------------|
| Methyl Ethyl Ketone (MEK)                            |                 |
| Cyclohexanone/MEK 90/10                              |                 |

| Solvent bonding to CYROLITE® and other medical grade resins using Cyclohexanone/MEK 90/10 |                 |    |      |
|-------------------------------------------------------------------------------------------|-----------------|----|------|
| Tubing materials                                                                          | CYROLITE® GS-90 | PC | MABS |
| Silicone                                                                                  |                 |    |      |
| Polybutene                                                                                |                 |    |      |
| non-DEHP Poly Vinyl Chloride (PVC)                                                        |                 |    |      |
| *Thermoplastic Polyurethane (TPU)                                                         |                 |    |      |

\* Excellent bonding to TPU


| Legend                        |       |
|-------------------------------|-------|
| Lap shear pull force strength | Color |
| Poor bond                     |       |
| Good/Excellent bond           |       |

# Hansen Solubility parameter - Screening for compatibility - Substrate, tubing, solvent combinations



| *Hansen Solubility Parameter $\delta$ ((MJ/m <sup>3</sup> ) <sup>1/2</sup> ) |            |            |            |
|------------------------------------------------------------------------------|------------|------------|------------|
| Substrate Polymer                                                            | $\delta_D$ | $\delta_P$ | $\delta_H$ |
| Acrylics                                                                     | 18.1       | 8.4        | 3.1        |

| *Hansen Solubility Parameter $\delta$ ((MJ/m <sup>3</sup> ) <sup>1/2</sup> ) |            |            |            |
|------------------------------------------------------------------------------|------------|------------|------------|
| Tubing material                                                              | $\delta_D$ | $\delta_P$ | $\delta_H$ |
| PVC                                                                          | 18.7       | 10         | 3.1        |
| TPU                                                                          | 18.8       | 10         | 8.2        |
| SBC                                                                          | 17.4       | 1.12       | 2.64       |
| Polyamide                                                                    | 18.5       | 5.1        | 12.2       |
| PTFE                                                                         | 14         | 0          | 0          |



| *Hansen Solubility Parameter $\delta$ ((MJ/m <sup>3</sup> ) <sup>1/2</sup> ) |            |            |            |
|------------------------------------------------------------------------------|------------|------------|------------|
| Solvents                                                                     | $\delta_D$ | $\delta_P$ | $\delta_H$ |
| Cyclohexanone                                                                | 17.8       | 8.4        | 5.1        |
| Methyl ethyl ketone                                                          | 15.9       | 8.7        | 5.1        |

# Conclusions

- Solvent bonding is a versatile technique to adhere tubing to polymer materials for medical interconnector applications
- CYROLITE® GS-90 presents a BPA-free, and latex-free bonding solution to non-DEHP PVC tubing as well as alternatives to PVC tubing
- Hansen solubility parameter as a viable tool for determining polymer (substrate-tubing) to solvent compatibility (reference slide 21)
- Cyclohexanone/MEK and MEK are suitable choices for optimum solvent bonding of CYROLITE® GS-90 to non-DEHP PVC tubing

Table 6.2. Hansen Solubility Parameters of Selected Polymers. Retrieved from <https://app.knovel.com/hotlink/itble/rcid:kpHPE00024/id:kt00XREBS1/handbook-plasticizers/table-6-2-hansen-solubility> Handbook of Solvents, Volume 1 - Properties (2nd Edition)

# Q & A



# CYROLITE®

ADVANCED  
MEDICAL ACRYLICS

Roehm America LLC  
299 Jefferson Road  
Parsippany, NJ 07834

Röhm GmbH  
Deutsche-Telekom-Allee 9  
64295 Darmstadt  
Germany

Roehm Chemical (Shanghai) Co., Ltd.  
55 Chundong Rd. Xinzhuang Industry Park  
201108, Shanghai  
China

[www.roehm.com](http://www.roehm.com)  
[www.cyrolite.com](http://www.cyrolite.com)

## RÖHM

TRADITIONALLY  
INNOVATIVE



COMING SOON: January 2021  
**Secondary Operations Part 2: Ultrasonic Welding**

Visit CYROLITE on LinkedIn for updates and registration

# RÖHM

TRADITIONALLY  
**INNOVATIVE**